

miniRaman Spectrometer

The world's smallest
Raman Spectrometer



- Premium performance for quantitative analysis
- Excellent for material identification
- Accurate with high sensitivity
- Extended spectral range via dual laser
- Affordable

Miraspec

Software for PC (Windows 11)
and smartphone (Android)

Data acquisition	Laser power control Exposure time control Sensor gain control Number of repetitions control Spectral range control
Data preprocessing options	Spike correction (Whittaker-Hayes, moving window) Spectral smoothing (Whittaker, asymmetric least squares, Savitzky-Golay) Baseline correction (rolling circle, rubberband, least squares, asymmetrically reweighted penalized) Spectrum normalization (Z-score, mean, Mean centre, MinMax, Unit Norm L1, Unit Norm L2) Intensity normalization* Laser mode-hop correction* Spectral super resolution*
Data Exploration	Principal component analysis Non-negative matrix factorization (SIMPLISMA-NNLS, MCR-ALS)
Qualitative analysis	Material identification (Pearson correlation, square Euclidean cosine, square first difference Euclidean cosine) Compatible libraries (>20,000 spectra)** Creation of spectral libraries
Quantitative analysis	Principal component regression Raman peak height/peak area calibration Partial Least Squares (PLS) calibration
Classification	Random Forest, Linear SVM, AdaBoost, Decision Tree, Hoefding Tree, Naïve Bayes, Perceptron, Softmax Regression

*patented feature

**various library options available upon request

Spectrometer system controlled by smartphone or PC via Bluetooth or USB-C cable.



Accessories

- Short/ Middle working distance probe, f=15 mm
- Long working distance probe, f=30 mm
- Contact probe
- Disposable silicone heads
- Sample/Vial holders
- Axial focusing accessories
- Light protection sample cover
- Power bank
- Belt Holder
- Docking station

Applications

For researchers, industries, universities and general consumers



Life Sciences & Health:

Biosciences, pharmaceuticals, skin diagnostics, cosmetics



Materials & Nanotechnology:

Polymers, nano-materials, semiconductors, surface enhanced Raman scattering (SERS)



Industrial & Chemical Analysis:

Chemicals, geology, forensics




Quality Control & Authentication:

Alcohol quality, counterfeit product detection



Unique Technology



The **miniRaman Spectrometer** is an innovative, miniaturized Raman spectrometer equipped with a patented, integrated **reference channel**. This innovative technology ensures **automatic calibration** by correcting for laser wavelength drift. The system continuously adjusts the Raman shift and intensity of the sample spectrum, delivering fast, accurate, and reliable results. These advanced features make the miniRaman spectrometer an ideal solution for **material identification** and precise quantitative measurements.

Specification

Feature versus model	miniRaman Standard	miniRaman SERS	miniRaman Dual	miniRaman Full Range	miniRaman Full Range SERS
Laser wavelength	785 nm		660 nm and 785 nm	660 nm	660 nm
Power on Sample*	5-50 mW	1-15 mW	5-45 mW (660) 5-50 mW (785)	10-80 mW	1-15 mW
Spectral Range	270-3100 cm^{-1}		2700-3700 cm^{-1} (660) 300-2500 cm^{-1} (785)	270-3700 cm^{-1}	270-3700 cm^{-1}
Spectral Resolution**	7, 9, 12 cm^{-1}		11, 14, 18 cm^{-1} (660) 7, 9, 12 cm^{-1} (785)	11, 14, 17 cm^{-1}	11, 14, 17 cm^{-1}
Signal-to-noise ratio***	500:1	100:1	600:1 (660) 500:1 (785)	500:1	100:1

* Actual laser power range can differ $\pm 5\%$ per device. Please contact us if you need specific laser power range values. Lower limit requires exposure time of 2 seconds

** Slit size dependent; slit size can be customized (20, 35, 50 μm slits)

*** Determined as peak signal-to-noise ratio of polystyrene spectrum at maximal laser power, integration time 0.3s, number of repetitions 10

Extra Features

- Small size (105 x 50 x 45 mm)
- Lightweight (360 g)
- One button operation
- Low power consumption (Up to 4 hr)
- Bluetooth & USB connectivity

Customisation

- Operator body - IP67 and Drop resistant body design



Contact Details

Lightnovo ApS

Blokken 11, 1.

3460 Birkerød

Denmark (DK)

+45 71 37 04 10

info@lightnovo.com

CVR: 40979603


About Lightnovo

A spin-off from the Technical University of Denmark, Lightnovo was founded in 2019 by an enthusiastic team united by the goal of revolutionizing the field of Raman spectroscopy through innovative, high-performance solutions. Our mission is to develop and commercialize “Raman for all: democratize the power of high-end Raman spectroscopy for the benefit of mankind”.

We aim to provide premium performance Raman spectrometers and microscopes with the world's smallest form factors without compromising the performance. With this innovation, Lightnovo addresses the need for portable, reliable field instruments at an affordable price.

© 2026 Lightnovo Aps. All rights reserved

Scan to learn more about this and other products



Please consult your local sales representative for details