

miniRaman microscope



User Manual

Version 001

For Technical Support, please contact:

Lightnovo ApS

Blokken 11, 1.

3460 Birkerød

Denmark

Phone: +45 71465729

<https://lightnovo.com/>

service@lightnovo.com

For International Support, please contact your local distributor.

Lightnovo ApS provides this document to its customers with a product purchase to use in the product operation. This document is copyright protected and any reproduction of the whole or any part of this document is strictly prohibited, except with the written authorization of Lightnovo ApS.

All technical information in this document is for reference purposes only. The contents of this document are subject to change without notice.

Lightnovo ApS makes no representations that this document is complete, accurate or error-free and assumes no responsibility and will not be liable for any errors, omissions, damage or loss that might result from any use of this document, even if the information in the document is followed properly.

For Research Use Only. This instrument or accessory is not a medical device and is not intended to be used for the prevention, diagnosis, treatment or cure of disease.

© 2026 Lightnovo Aps. All rights reserved.

Contents

Contents	3
Overview	4
Instrument description	5
General specifications	6
Configuration options	6
Application of use	7
Safety information	8
Qualified Personnel	8
Warning labels	8
Warning label positions	9
Safety	9
General instrument information.....	10
Safety Instructions	10
Installation	12
Scope of Delivery	12
Standard Components	12
Inspecting the Packaging	12
Transportation	12
Site requirements	13
Accessories	14
Software installation	16
Operation	17
Attaching & removing probes.....	17
Removing and replacing the spectrometer	19
LED guide.....	21
Hardware connection.....	22
Technical drawings	23
Service	25
Recommended care	25
Addresses	25

Overview

The miniRaman Microscope is an ultra-compact, high-performance confocal Raman microscopy system designed to deliver reliable, high-resolution chemical analysis in a flexible and accessible format. At the core of the miniRaman microscope is a patented integrated reference channel, which continuously corrects for laser wavelength drift and enables automatic calibration. This unique technology ensures excellent wavenumber accuracy, stable Raman intensity calibration, and exceptionally high optical throughput, resulting in accurate, reproducible, and comparable Raman measurements over time. The system supports high-resolution Raman mapping, allowing detailed spatial characterization of heterogeneous samples, while its dual operation modes enable seamless use of the spectrometer both integrated within the microscope and as a standalone Raman unit. The modular architecture provides flexibility in configuration and future upgrades, adapting easily to different applications and laser wavelengths.

Combining compactness, performance, and affordability, the miniRaman microscope delivers confocal Raman microscopy without compromise, making advanced Raman analysis accessible to a wider range of laboratories and users.

Instrument description

The miniRaman microscope is a compact portable Raman microscope. Physical dimensions of approximately 175 mm x 180 mm x 410 mm, weighing approximately 4 kg. The instrument has 2 USB C ports on the bottom of the back, one for charging and the other for communication.



Figure 1. Front, back, top and bottom view of miniRaman microscope

General specifications

Table 1. Typical values for the specifications of the miniRaman microscope.

Feature versus model	miniRaman microscopy platform motorized	miniRaman microscopy platform manual	miniRaman microscopy platform focusing stage
Compatible spectrometers	miniRaman spectrometers miniRaman Pro spectrometers		
Lateral resolution*	800 - 900 nm		
Axial resolution or confocality*	2.5 – 3 μ m		
White light microscopy	Reflected with simultaneous visualization of laser spot and Raman acquisition		Not included
Microscopy configuration	up-right and inverted		
Mapping travel range in XYZ	30 x 30 x 30 mm	10 x 10 x 13 mm	
Minimum step size	500 nm	5 μ m	
Physical dimensions	175x180x410 mm (LxWxH)		
Weight	11.5 kg	10.5 kg	9.5 kg

* Resolution depends on the laser wavelength and spectrometer used. Values represent the use of miniRaman spectrometers, all wavelengths with microscope objective NA = 0.95, 100x magnification

Configuration options

Operator

The Operator configuration adds water and dust proofing, as well as drop resistance to the instrument.

Configuration matrix

The following table explains which configuration options are compatible with which laser versions on the miniRaman microscope.

Table 2. Configuration options for miniRaman microscope.

	mRm motorised	mRm manual	mRm focusing stage
Operator	Yes	Yes	Yes



Note

All configurations' options for the miniRaman spectrometers see in miniRaman manual.

Application of use

The miniRaman microscope is a general purpose compact Raman microscopy analysis instrument for qualitative and semi-quantitative analysis of liquids and solids. It is intended for laboratory, industrial and field environments. Raman spectroscopy is non-destructive and can be performed without contact with the sample. This versatile and modular instrument can be used for identification, verification, and screening of a wide variety of samples, using library searching. Semi-quantitative and quantitative measurements can be performed using chemometric analysis on Lightnovo's computer software. The instrument also has the capability to perform mapping measurements of stationary samples.

Note



In Raman spectroscopy, sample fluorescence can yield a much more intense signal than the Raman scattering signal of the sample, masking any Raman bands in the spectrum. Therefore, Raman spectroscopy is normally not a suitable analysis technique for fluorescent samples.

Safety information

Read the following safety information carefully before operating the Lightnovo instrument and keep this manual for future reference available at any time. Always follow the instructions described in this manual to ensure safety operation and to avoid property damage. Improper use or failure to follow these safety instructions can result in serious injuries and/or property damage. Any non-observance of the precautions will infringe the intended use (i.e. performing measurements by Raman spectroscopy) of the instrument. In this case Lightnovo ApS will not assume any liability. It is the operator's duty to plan and implement all necessary safety measures and to supervise their observance. Moreover, the operator must ensure that the instrument is in proper functioning condition. A safe and faultless operation can only be guaranteed if the Lightnovo instrument is transported, stored, installed, operated and maintained properly according to the procedures described in this manual. Never remove or deactivate any safety systems during the instrument's operation. Ensure that objects and/or material not required for the measurement is out of the instrument operating area.

Qualified Personnel

Primary installation maintenance and repair works not described in this manual should only be performed by Lightnovo approved service personnel. Only authorized operating personnel that have been briefed about the instrument operation and all relevant safety aspects should operate and maintain (i.e. only maintenance works that are described in this manual) the instrument. It is the duty of the instrument owner to ensure users have proper training for use of a laser instrument and relevant risk assessments are in place. All repairs, adjustments and alignments on any instrument component must be performed in accordance with the safety regulations and standards applied in the country in which the instrument is installed.

Warning labels

This warning symbol indicates the existence of laser radiation.



Class 3B lasers are hazardous for eye exposure. They can heat skin and materials but are not considered a burn hazard. For visible-light lasers, Class 3B lasers' output power is between 5 and 499 milliwatts. Class 3B lasers are normally hazardous under direct beam viewing conditions, but are normally safe when viewing diffuse reflections.

Warning label positions

A warning label can be found on the rear of the miniRaman microscope. A warning symbol can also be found next to the port on the top of the instrument where the laser radiation exits the instrument.



Figure 2. Warning label for miniRaman microscope

Safety

The following chapters describe all relevant safety aspects of the instrument operation. Depending on the degree of hazard the safety instructions are classified as follows:



Danger

indicates that death, severe personal injury or substantial property damage **WILL** result if proper precautions are not taken.



Warning

indicates that death, severe personal injury or substantial property damage **CAN** result if proper precautions are not taken.



Caution

indicates that minor personal injury or property damage **CAN** result if proper precautions are not taken. Important draws your attention to a particularly important information.



Note

draws your attention to useful information on the product, e.g. product operation or to a special part of the manual.

The safety instructions Danger, Warning and Caution are marked by the corresponding warning labels.

General instrument information

The miniRaman microscope uses a laser (660 and/or 785nm wavelength). The laser diodes used in the spectrometer emit continuous visible and partially visible laser radiation in the near infrared region. According to the standard EN 60825-1:2014, miniRaman spectrometer is Class 3B laser product. Therefore, it is normally hazardous under direct beam viewing conditions but is normally safe when viewing diffuse reflections.

The door on the front of the instrument can have an enabled interlock, making the instrument a Class 1 laser product.

Safety Instructions

In addition to the safety instructions given below, one should also comply with all local regulations concerning laser safety.

The analysis system is specified as a laser class 3B product if the interlock is not enabled, i.e. it is considered safe if handled carefully, with restricted beam viewing.

It is advised that laser safety goggles that meet the specifications for the laser your system uses are worn by the user and others within the vicinity. Laser glasses should have an LB rating above 3 for code D.

Warning



Avoid eye and skin exposure to direct or scattered laser radiation under all circumstances! Failure to do so can cause permanent and irreversible eye damage and/or serious skin injuries!

Below is a table of the nominal ocular hazard distance (NOHD) values for the notable configurations at maximum laser power. Be aware these are calculated assuming direct line of sight to the objective, which is not possible in this microscope. However, these values could be useful guidance if the samples reflect the laser out of the sample compartment.

Table 3. NOHD values for the mRs microscope standard configuration with different laser configurations (columns) using different probes (rows) at max laser power. Note that the actual maximum laser power of the instrument may differ from these values.

mRs wavelength	660 nm	785 nm
Maximum power	80 mW	50 mW
10x objective	0.26 m	0.42 m
20x objective	0.09 m	0.13 m
50x objective	0.25 m	0.4 m
100x objective	0.05 m	0.06 m
100x NIR objective	0.06 m	0.07 m

Installation

Scope of Delivery

Standard Components

- miniRaman spectrometer (including user manual and quality test report)
- miniRaman microscope (including user manual and quality test report)
- Accessories (includes spares, adaptors, objective lenses, sample preparation tools etc.) typically:
 - Objective adaptor
 - Microscope objective
 - Goggles for laser protection
 - Accessories for independent spectrometer use
 - Accessories for sample preparation
- USB-C data cable
- USB-C power adaptor

Inspecting the Packaging

After having received miniRaman microscope, inspect the packaging for damages. If there are any signs of damage, contact your local shipping representative before opening the shipping box.

Warning



Do not put miniRaman spectrometer and microscope into operation if there are signs of damage. Failure to do so may result in severe personal injuries and/or property damage.

Transportation

When transporting the microscope, we recommend you use the Lightnovo backcarry on bag to avoid damages.

When transporting the miniRaman microscope by plane it is recommended to use the Lightnovo transport safety trolley.



Figure 3. Transportation case for miniRaman microscope.

Site requirements

Space Requirements

miniRaman microscope requires a space of 20 cm² benchtop space. (For the exact instrument dimensions refer to Specification.) At the rear instrument side, take a clearance of at least 3cm into account.

Environmental Requirements

To ensure optimal instrument performance and long-term reliability the following environmental conditions are essential:

- Ambient temperature range: 18 - 35°C (64 - 95°F)
- Humidity (non-condensing): ≤ 70% (relative humidity)



Note

miniRaman spectrometer is an instrument of protection class I (electrical safety).

Accessories

One will receive in the package a set of accessories. Depending on the order each accessory is recommended to be used for a particular type of measurement and test.

The probes for contact and distance measurements should be oriented by the upper side as shown in Table 4.

Table 4. Accessories for Lightnovo microscope

Accessory	Description	Type of measurement and test
	10x objective	ZEISS Olympus
		
	20x objective	Olympus
	50x objective	Olympus
	100x objective	ZEISS

	100x NIR objective	ZEISS
	RMS Thread adaptor for Olympus objective	To attach the objective to the magnetic adaptor
	M25 adaptor for Nikon objectives	To attach the objective to the magnetic adaptor
	M27 adaptor for Zeiss objective	To attach the objective to the magnetic adaptor
	Magnetic adaptor	To attach the objectives to the microscope
	M27/0.75 thread adaptor for LN Spectrometers for Zeiss objectives M25/0.75 thread adaptor for LN Spectrometers for Nikon objectives (+ extra metal adaptors)	To attach the objectives to the microscope

Software installation

1. Download Miraspec software package from the Lightnovo website:
<https://lightnovo.com/miraspec/>
2. Open the manual inside for detailed installation and instrument connection instructions.

Operation

Attaching & removing probes

To attach an objective, check the threading of the objective.



If the objective is M27 threaded, go to the next step.



If the objective is M25 or RMS, then first attach the adaptor ring to the objective.



Screw the objective (with adaptor if needed) onto the magnetic adaptor (M27).



Place the objective with the magnetic adaptor at the top into the slot in the microscope. The magnetic adaptor may need to be rotated to lock into place.



Removing and replacing the spectrometer

Press the button at the top on the front of the microscope, this will open the top hatch instrument.



Open the top hatch all the way.



Push the bar at the top back to disengage the lock position.



Pull the bar until it can relax backwards, and the spectrometer is visible.



Disconnect the USB-C connect to the instrument.



Open the two hinged holders of the instrument.



Grab the spectrometer from the sides and pull it free instrument.



LED guide

Blue

The microscope is connected.



Blue/ Green

The system is connected to computer, the laser is on and the system is ready to operate.



Blue/Green/Red

The system is actively mapping.



Blue / Red

The system is turning off.



Hardware connection


To connect the microscope to a computer first plug the USB C cable into the port on the bottom of the rear of the microscope.



Plug the USB end of the cable into the computer.



Open the Lightnovo software on the computer. Press the connect button on the software to begin communication with the instrument. (see software manual for details).



Technical drawings




Figure 4. miniRaman microscope drawing with dimensions

Service

Recommended care

Switching the system ON/OFF

When the analysis system is not used for a longer period of time, it is highly recommended to switch off the Raman excitation laser. This action will prolong the service lifetime of the laser.

Warning



Do not work with miniRaman spectrometer at laser power that is higher than specified for particular model. This could lead to the decreased life time of the laser diode or laser damage.

Addresses

Lightnovo ApS

Blokken 11, 1.

3460 Birkerød

Denmark

Phone: +45 71465729

<https://lightnovo.com/>

service@lightnovo.com



lightnovo.com