



# Precise Focusing Stage

Accessory could be combined with all models of miniRaman spectrometer

- Provides protection from ambient light
- Provides a safe operation since user will not see the laser irradiation
- Provides precise focusing of the laser beam on a sample

## Precise Focusing Stage XYZ



- lateral travel range: 25 mm along the X and Y axes
- axial travel range: 25 mm
- precise focusing of the laser beam on a sample in axial and lateral dimensions with accuracy 10  $\mu\text{m}$
- magnetic fixation of the stage cover in up and down positions

## Precise Focusing Stage Z



- axial travel range: 25 mm
- precise focusing of the laser beam on a sample in axial dimension with accuracy 10  $\mu\text{m}$
- magnetic fixation of the stage cover in up and down positions



**Lightnovo ApS**  
Blokken 15, 1. tv.  
3460 Birkerød  
Denmark (DK)  
CVR: 40979603

+45 71 37 04 10  
[info@lightnovo.com](mailto:info@lightnovo.com)  
<https://lightnovo.com>

