

# Precise Focusing Stage



**Accessory can be combined with all models of miniRaman spectrometer**

- Provides protection from ambient light
- Provides a safe operation since user will not see the laser irradiation
- Provides precise focusing of the laser beam on a sample

## Precise Focusing Stage XYZ



- lateral travel range: 25 mm along the X and Y axes
- axial travel range: 25 mm
- precise focusing of the laser beam on a sample in axial and lateral dimensions with accuracy 10  $\mu$ m
- magnetic fixation of the stage cover in up and down positions
- weight: 3,2 kg
- dimensions: 215 mm x 175 mm x 320,5 mm

## Precise Focusing Stage Z



- axial travel range: 25 mm
- precise focusing of the laser beam on a sample in axial dimension with accuracy 10  $\mu$ m
- magnetic fixation of the stage cover in up and down positions
- weight: 2,5 kg
- dimensions: 215 mm x 175 mm x 320,5 mm



**Lightnovo ApS**  
Blokken 15, 1. tv.  
3460 Birkerød  
Denmark (DK)  
CVR: 40979603

+45 71 37 04 10  
[info@lightnovo.com](mailto:info@lightnovo.com)  
<https://lightnovo.com>

